
North American Distribution:  
**AMD Technologies, Inc.**  
V/F: 1.310.471.8900  
giesberg@amdtechnologies.com  
218 Bronwood Avenue, Los Angeles, CA 90049-3104 USA

# pehamed

## DINo 152,7

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The DINo 152,7 phantom can be used to perform both the  
acceptance and the constancy tests of conventional  
mammography exposure units.

The phantom fulfills all of the requirements of the DIN V 6868, Part 152  
and the DIN 6868, Part 7 standards.

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## DINo 152,7


Using the DINo 152,7 phantom, the following parameters of mammography exposure units and the corresponding film processors can be determined.

The upper detail plate is fixed with pins to permit easy interchangeability with future detail plates, e.g. for digital mammography.

The MAMMOSOFT Module 152.7 software offers convenient administration, display and printing of the stored data.

- 1 optical density at the reference point
- 2 spatial resolution: lead bar patterns perpendicular and parallel to the chest wall side (8, 10, 13, 16 lp/mm).
- 3 contrast: 3 different absorption levels (cut-out, 3 mm deep step in the structure plate, 3.5 mm acrylic glass + 2.5 mm polytetrafluorethylene (PTFE, Teflon)).
- 4 threshold contrast: 7 round drill holes,  $\varnothing$  5.5 mm, depths: 0.1 mm, 0.15 mm, 0.2 mm, 0.25 mm, 0.3 mm, 0.35 mm, 0.4 mm.
- 5 boundary of the chest wall side image: 4 rows of spheres, each consisting of 5 adjacent steel spheres,  $\varnothing$  2 mm.
- 6 film processing (Dmin, D step 21, LE, LK, etc.), 0.1 mm lead foil prevents film exposure.

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The phantom system is made up of 3 parts: a 40 mm and a 20 mm thick absorption plate, and a 6 mm thick structure plate. The phantom has the following overall dimensions:

180 x 240 x 46 mm

Pehamed, Geräte GmbH, Mühlstrasse 38, D - 65843 Sulzbach/Ts. (Germany / Europe) Tel.: +49 6196 500 40 30  
Fax: +49 6196 500 40 50 E-mail: [pehamed@t-online.de](mailto:pehamed@t-online.de) Internet: <http://www.pehamed.com>